Dale® ACE Connector® with Legacy Connection

Infuse nutrients, medications, fluids, irrigate, check gastric residuals, control flow and suction

APPLICATION INSTRUCTIONS

Product Number
475



ON/OFF HANDLE

controls fluid flow.

"ON" for feeding or drainage.

"OFF" to irrigate tube, check gastric residuals, administer medication, or ambulate patient.

PROXIMAL END

(Feeding Bag End)

furthest from patient.

- connects to enteral pump and gravity delivery sets.
- connects to bolus feeding syringes.
- connects to wall or portable suction.

SYRINGE MED PORT SEAL

accepts most 50-60cc catheter-tip syringes for irrigation, medication delivery and checking gastric residuals.

"Push & Twist" Syringe or Feeding Administration Set stepped connector into syringe med port.

REMOVABLE WHITE FLEXIBLE ADAPTER

When connecting to suction or drainage tubing

DISTAL END

(Patient End)

closest to patient's stomach.

- universally connects to most enteral tubes.

CATHETER CONNECTOR CONNECTS TO:

- Most small bore feeding tubes
- Salem Sump™ tubes
- Levin tubes
- PEG and PEJ tubes*
- Balloon replacement PEG & PEJ tubes*
- Most G-J tubes*

*18FR and higher

THE ACE ADVANTAGE

- Supports a closed enteral system
- Protects clinician from potentially infectious splash-back
- Protects patient from outside contamination
- No soiled bed linens or loss of formula due to open adapters
- Provides auto-seal syringe med port seal for accepting medication, irrigation, or gastric residuals

Note: the distal end is always on the side closest to the pt. This is true from a clinician's perspective. It is exactly opposite from a patient's point of view.

Dale Medical Products

100% EMPLOYEE OWNED AND OPERATED

800-343-3980 • www.dalemed.com

SCAN FOR MORE INFO

Dale & ACE Connector are registered trademarks of Dale Medical Products, Inc. in the USA, EU, and UK.
Salem Sump is a trademark of Cardinal Health, Inc.
©2024 Dale Medical Products, Inc. All rights reserved.