Dale
Medical Products

Product Catalog

Established in 1961, Dale is recognized as a trusted manufacturer of specialty medical devices that provide high quality, reliable, cost-effective solutions to enhance patient care.





Abdominal Binders

Product #: 410, 411, 418, 525, 810, 811, 818, 820, 920



ACE Connectors®

Product #: 475, 485




Bendable ArmBoards

Product #: 650, 651, 652, 653



Endotracheal Tube Holders

Product #: 270, 280




Hold-n-Place® Foley Catheter Holders

Product #: 316, 330




IV-ARMOR®

Product #: 675



IV Line Manager

Product #: 630



Multipurpose Tube Holder

Product #: 150



Nasal Dressing Holder

Product #: 600



NasoGastric Tube Holder

Product #: 160, 161, 163



Post-Surgical Bras

Product #: 701, 702, 703, 704, 705, 706, 712, 713, 714, 715, 716



Scrape-n-Suction® Tongue Cleaner

Product #: 220



Tracheostomy Tube Holders

Product #: 240, 241, 242, 244



Transducer Holders

Product #: 620, 622

Quality devices you can count on


Abdominal Binder with EasyGrip™ Strip

Dale Abdominal Binders are the only binders clinically shown to expedite the return of pulmonary function after open abdominal surgery while encouraging patients to ambulate sooner. Available in sizes up to 94” long and in widths up to 15”.

- Exclusive stretch material provides full, all-around compression.
- Hook & loop closure secures at any point along the binder circumference.
- Promotes mobility and lessens the chance of sedentary complications.
- Available in bariatric sizes.
- EasyGrip™ Strip holds up to four 100cc drainage bulbs for easier monitoring and changing.







Product #	Description	Quantity
525	**Abdominal Binder, can also be used a rib belt** 6” Wide, Stretches to fit 72cm-127cm (28”-50”)	1/box
410	**Abdominal Binder** 9” Wide, Stretches to fit 76cm-114cm (30”-45”)	1/box
411	**Abdominal Binder** 9” Wide, Stretches to fit 117cm-157cm (46”-62”)	1/box
418	**Abdominal Binder** 9” Wide, Stretches to fit 152cm-191cm (60”-75”)	1/box
810	**Abdominal Binder** 12” Wide, Stretches to fit 76cm-114cm (30”-45”)	1/box
811	**Abdominal Binder** 12” Wide, Stretches to fit 117cm-157cm (46”-62”)	1/box
818	**Abdominal Binder** 12” Wide, Stretches to fit 152cm-191cm (60”-75”)	1/box
820	**Abdominal Binder** 12” Wide, Stretches to fit 183cm-239cm (72”-94”)	1/box
920	**Abdominal Binder** 15” Wide, Stretches to fit 183cm-239cm (72”-94”)	1/box

Patents dalemed.com/patents

ACE* Connectors®

*Access Controller for Enteral

The Dale ACE Connector (Access Controller for Enteral) is the simple and secure way to connect tubes for in-line feeding, flushing, irrigation, sampling gastric residuals and medication delivery.

- 100% closed system reduces the risk of potentially infectious splash back and cross-contamination.
- On-off handle controls fluid flow—turn on for feeding or off for tube irrigation, medication administration, gastric residuals, or patient ambulation.
- Fully integrated auto syringe seal med port eliminates splash back and can be used for medication delivery, irrigation, and sampling of gastric residuals. ENFit compatible ACE accepts most ENFit size syringes from .5 mL to 100 mL. The Legacy ACE will accept most 50 and 60 mL catheter tipped syringes for the same functions.
- Connects to most NG, feeding tubes, Salem Sumps™, Suction tubes, and other enteral tubes, including most PEG, PEJ, and G-J tubes.







Original 475 ACE Legacy Connector

ENFit-80369-3 compliant 485 ACE Connector

MR Conditional

Product #	Description	Quantity
475	**ACE Connector®** Fits most Legacy NG enteral feeding tubes, Legacy stomach tubes, and 50-60mL catheter tip syringes	50/box
485	**ACE Connector® with ENFit® Connections** Fits most ENFit NG enteral feeding tubes, Enteral ENFit stomach tubes, and .5-100mL ENFit tip syringes	50/box

Not available for purchase in Canada & EU

Patents dalemed.com/patents

Bendable ArmBoards

The Dale Bendable ArmBoard allows the patient's hand and wrist to maintain their natural position while securely and comfortably preventing catheter movement.



- Bendable inner core can be custom shaped to desired position – Enables clinician to hyperextend patient's wrist for pronounced radial artery exposure, trouble-free insertion and line integrity.
- Adjustable hook & loop closures provide comfortable fit.
- Elastic straps with adjustable hook & loop closures provide a comfortable fit.
- Four convenient sizes.
- MR Conditional.





Product #	Description	Quantity
650	**Bendable ArmBoard, Large** 23cm x 9cm (9" x 3.5")	10/box
651	**Bendable ArmBoard, Medium** 13.7cm x 6.4cm (5.38" x 2.5")	10/box
652	**Bendable ArmBoard, Small** 11.4cm x 4.4cm (4.5" x 1.75")	10/box
653*	**Bendable ArmBoard, X-Small** 10.8cm x 2.5cm (4.25" x 1")	10/box

*Detachable Straps.

Endotracheal Tube Holders

Dale Endotracheal Tube Holders provide reliable stabilization of endotracheal tubes(sizes 7.0-10.0mm) following intubation.

- Low profile, cushioned neck band is soft and comfortable on the face and neck.
- Skin-friendly adhesive can be left in place up to 3 days.
- Tubing channel is made of lightweight, flexible material to securely hold the endotracheal tube.
- Hook & loop closures provide a secure fit.
- Easy repositioning reduces risk of skin and lip ulcerations.
  *Tube Track on BreezeLock
- Easy access for oral care.
- Can be used in supine or prone position.





280 BreezeLock



270 Stabilock

Product #	Description	Quantity
270	**Stabilock** Fits endotracheal tube sizes 7.0mm-10.0mm	10/box
280	**BreezeLock® with Tube Track** Fits endotracheal tube sizes 7.0mm-10.0mm	10/box

Patents dalemed.com/patents


Hold-n-Place® Foley Catheter Holders

Dale Hold-n-Place Foley Catheter Holders stabilize the indwelling catheter, reducing the likelihood of traumatic removal, without restricting patient's movement or irritating the skin. All Dale Hold-n-Place Foley Catheter Holders feature unique locking systems that securely anchor the catheter in place, reducing the risk of urethral erosion, bladder spasms and trauma.

- The leg band, made entirely of soft stretch material, distributes compression evenly around the thigh to minimize circumferential pressure. Ideal for long-term catheterization.
- The multipurpose tube holder, with a skin-friendly adhesive base, is ideal for short-term catheterized patients.
- The locking devices speed up application and simplify adjustment.
- The holders can be used with equal effectiveness on all Foley catheters including Teflon, silver and silicone coated catheters.
- Leg band/waist band ideal for bariatric patients.



Hold-n-Place Leg Band



Hold-n-Place Waist Band



Hold-n-Place Multipurpose Tube Holder



Product #	Description	Quantity
316	**Hold-n-Place® Leg Band** Fits up to 76cm(30")	10/box
330	**Hold-n-Place® Leg Band/Waist Band** Fits up to 142cm(56")	10/box
150	**Hold-n-Place® Multipurpose Tube Holder** Fits catheters and tubes 12fr-40fr	50/box

Patents dalemed.com/patents


Hold-n-Place® Multipurpose Tube Holder


The Dale Hold-n-Place Multipurpose Tube Holder has a skin friendly, breathable, adhesive base that secures a wide variety of catheters and tubes from 12fr up to 40fr, including: Jackson-Pratt, PEG-tubes, suprapubic tubes, IV lines, chest tubes, Foley catheters, drainage tubing, and multi-port catheters.

- Compatible with a wide range of catheters and tubes (12Fr–40Fr), including Teflon™, silver, and silver-coated options for standardized use.
- Skin-friendly, breathable adhesive base allows secure wear for up to 7 days.
- Easy-to-apply base with unique closure for faster application and simple adjustments.
- Offers choice of adhesive or non-adhesive tube retention.
- DEHP-free and not made with natural rubber latex.









Product #	Description	Quantity
150	**Hold-n-Place® Multipurpose Tube Holder** Fits catheters and tubes 12fr-40fr	50/box

Patents dalemed.com/patents

IV-ARMOR®

Dale's IV-ARMOR minimizes downstream occlusions while enabling patient movement. The flexible protective overlay is the ideal way to maximize IV patency while protecting the site from patient tampering.

- Reduces the introduction of infection associated with IV reinsertion.
- Reduces downstream occlusion alarms.
- Reduces risk of catheter kinking in a bend such as a wrist or elbow.
- Reduces risk of catheter movement in an area of bony prominence, such as back of hand or forearm.
- Skin friendly adhesive.










Product #	Description	Quantity
675	IV-ARMOR® One size fits most	10/box


Patents dalemed.com/patents

IV Line Manager

The IV Line Manager organizes up to six (6) IV tubing lines into dedicated retention channels. By keeping lines separated and secure, the IV Line Manager prevents tangles which reduces the risk of infusion confusion and medication errors, and enables accurate, timely line tracing—a critical factor in clinical care. Per AAMI guidelines, IV lines should be organized by separating IV infusions and minimizing tangles.

- Holds and organizes up to six IV lines securely in place for streamlined workflow.
- Up to three IV Managers can be positioned—at proximal, intermediate, and distal points—for maximum control and organization.
- Prevents kinks, dislodgements, and tangling of lines with a clean, organized structure.
- Facilitates precise and efficient tracing of IV infusion lines from patient to source, supporting safe medication and fluid administration.
- Eliminates time wasted on detangling lines.
- Reduces risk of HAI by keeping lines organized and reducing frequent handling of lines which increases likelihood of contamination.





Product #	Description	Quantity
630	**IV Line Manager** organizes up to 6 IV lines	10/box

Patents dalemed.com/patents

Nasal Dressing Holder

The Dale Nasal Dressing Holder replaces adhesive tape to hold nasal dressings in place following rhinoplasty, septoplasty, sinus surgery, facial trauma and nose bleeds. Clear plastic dressing holder is secured by ear loops. Hook and loop fastener tabs can be adjusted to fit any patient comfortably.

- Soft moisture repellent material adds to comfort.
- Plastic center section holds dressings in place without absorbing secretions.
- Easy to remove and replace, making dressing changes simple and painless.





Product #	Description	Quantity
600	**Nasal Dressing Holder** One size fits most	10/box

Nasogastric Tube Holder

The Dale Nasogastric Tube Holder is ideal for holding most nasogastric tubes, including nasoenteric, Salem Sump, and Levin stomach tubes.

- Skin-friendly adhesive stays in place for up to 3 days, yet is easy to remove.
- Breathable material.
- Unique dual interlocking tabs(160 & 161) or single tab(163) spiral around the nasogastric tube for superior tube security.
- Non-adherent plastic tips make tube removal safer and easier, especially when wearing gloves.
- Single patient use only.
- No scissors required.
- May also be used to secure most manometry catheter and pH impedance probes used for esophageal manometry and pH impedance testing.





160 NasoGastric Tube Holder



163 UltraGrip NasoGastric Tube Holder

Product #	Description	Quantity
160	**NasoGastric Tube Holder** Large/Adult, Dual Tabs	50/box
161	**NasoGastric Tube Holder** Small/Young Adult/Pediatric/Geriatric, Dual Tabs	50/box
163	**UltraGrip® NasoGastric Tube Holder** Large/Adult, Single Tab	50/box

Patents dalemed.com/patents

Post-Surgical Bras

The Dale Post-Surgical Bras provides the support and compression patients need immediately following breast and thoracic surgery. Soft, breathable fabric stretches to accommodate edema without binding or irritating sensitive tissue.

- Front hook & loop closure facilitates inspection of the surgical site and dressing changes.
- Adjustable shoulder straps allow for ease of application and improved comfort
- Can be used to hold surgical dressings in place without tape.
- Indicated for use after a wide range of surgical procedures: Reduction, reconstruction, lift, mastectomy, lumpectomy, thoracic, biopsy, and fine needle aspiration.
- Ideal for post radiation treatment.
- Draingage bulb holding tabs.*
  *Compression Bras only



Post-Surgical Bra Light Compression	Product #	Sizes	Quantity
	701	**Small** fits A-C 81cm - 86cm (32" - 34")	1/package
	702	**Medium** fits B-D 86cm - 92cm (34" - 36")	1/package
	703	**Large** fits B-D 92cm - 97cm (36" - 38")	1/package
	704	**X-Large** fits B-D 97cm - 112cm (38" - 44")	1/package
	705	**XX-Large** fits C-E 117cm - 137cm (46" - 54")	1/package
	706	**XXX-Large** fits E-I 137cm - 152cm (54" - 60")	1/package

Post-Sugical Compression Bra Moderate Compression	Product #	Sizes	Quantity
	712	**Medium** fits B-D 86cm - 92cm (34" - 36")	1/package
	713	**Large** fits B-D 92cm - 97cm (36" - 38")	1/package
	714	**X-Large** fits B-D 97cm - 112cm (38" - 44")	1/package
	715	**XX-Large** fits C-E 117cm - 137cm (46" - 54")	1/package
	716	**XXX-Large** fits E-I 137cm - 152cm (54" - 60")	1/package

Scrape-n-Suction® Tongue Cleaner

The Dale Scrape-n-Suction Tongue Cleaner scrapes the tongue with suction to remove bacteria and toxins in the tongue's biofilm during the oral care routine.

- Effectively removes thick, encrusted biofilm from the tongue.
- Tongue scraping is more effective at eliminating biofilm than tooth brushing alone.[1,2,3]
- May help reduce the risk of ventilator-associated pneumonia (VAP) and other respiratory infections, improving overall patient outcomes.[4,5]
- Designed to reach the back of the tongue with ease.
- Compact size with a uniquely curved front end for enhanced maneuverability.
- Fits comfortably in the oral cavity, even when occupied by an endotracheal tube (ET Tube).
- Easily attaches to suction for convenient removal of scraped biofilm.







Product #	Description	Quantity
220	**Scrape-n-Suction® Tongue Cleaner** Large/Adult	50/box

Patents dalemed.com/patents

[1]Kuo, Y. W., Yen, M., Fetzer, S., & Lee, J. D. (2013). Toothbrushing versus toothbrushing plus tongue cleaning in reducing halitosis and tongue coating: a systematic review and meta-analysis. Nursing research, 62(6), 422–429. https://doi.org/10.1097/NNR.0b013e3182a53b3a

[2]Santos, P. S., Mariano, M., Kallas, M. S., & Vilela, M. C. (2013). Impact of tongue biofilm removal on mechanically ventilated patients. Revista Brasileira de terapia intensiva, 25(1), 44–48. https://doi.org/10.1590/s0103-507x2013000100009

[3]Bordas A, McNab R, Staples AM, Bowman J, Kanapka J, Bosma MP. "Impact of different tongue cleaning methods on the bacterial load of the tongue dorsum." Archives of Oral Biology. 2008.

[4]Izumi M, Akifusa S. "Tongue cleaning in the elderly and its role in the respiratory and swallowing functions: Benefits and medical perspectives." Journal of Oral Rehabilitation. 2021 Dec;48(12):1395-1403.

[5]Takahama, A., Jr, de Sousa, V. I., Tanaka, E. E., Ono, E., Ito, F. A. N., Costa, P. P., Pedriali, M. B. B. P., de Lima, H. G., Fornazieri, M. A., Correia, L. S., Cardoso, L. T. Q., & de Maio Carrilho, C. M. D. (2021). Analysis of oral risk factors for ventilator-associated pneumonia in critically ill patients. Clinical oral investigations, 25(3), 1217–1222. https://doi.org/10.1007/s00784-020-03426-x

Tracheostomy Tube Holders

The Dale Tracheostomy Tube Holder is designed to provide secure positioning and minimize movement of the tracheostomy tube. Unnecessary movement of the trach tube can result in accidental dislodgement or displacement of the tube, tracheal fistula, tracheal stenosis, or airway granuloma.

- Eliminates the frustrations associated with twill ties and other holders while minimizing secondary complications.
- Moisture wicking neck band lining reduces the risk of skin breakdown.
- Fastener tabs quickly and securely attach to the tracheostomy tube holder.










Product #	Description	Quantity
240	**240 Blue®** 2.54cm(1") wide neckband, fits up to 49.5cm(19.5") neck	10/box
241	**241 PediPrints®** 2.54cm(1") wide neckband, fits up to 46cm(18") neck	10/box
242	**242 PediPrints®** 1.91cm(0.75") wide neckband, fits up to 23cm(9") neck	10/box
244	**244 Blue®** 3.81cm(1.5") wide neckband, fits up to 71cm(28") neck	10/box


Patents dalemed.com/patents

Transducer Holders

The Dale Transducer Holders offers a level and secure platform for up to three transducers. Stability ensures correct waveform readings and corresponding numerical display.

- Platform conforms to patient's arm.
- Soft, comfortable material with adjustable hook and loop closures.
- Secures up to three transducers.
- Holding tabs reduce the time associated with manual re-leveling and taping individual transducers to the patient's chest, arm, or mounting on the bedside pole.
- Increases organization and security of hemodynamic catheters during patient transportation and ambulation.





Product #	Description	Quantity
620	**Transducer Holder, fits up to three standard transducers** One size fits most	10/box
622	**Transducer Holder, fits up to two large transducers** One size fits most	10/box

Dale® Medical Products

508-695-9316

www.dalemed.com

©2026 Dale Medical Products, Inc. All rights reserved.
Dale, ACE Connector, Hold-n-Place, 240 Blue, 244 Blue, PediPrints, BreezeLock, IV-ARMOR, UltraGrip, & Scrape-n-Suction are registered trademarks of Dale Medical Products, Inc. in the USA, EU, and UK.
EasyGrip is a trademark of Dale Medical Products, Inc. in the USA.
Salem Sump is a trademark of Cardinal Health, Inc.
ENFit is a registered trademark of GEDSA used with their permission.
AD-001_REV_Sw

Scan or Visit dalemed.com to learn more